

OXYGEN ELECTROCHEMICAL PUMP

GEN' AIR



Based on the ionic conduction properties of zirconia, **GEN' AIR** makes it possible to create and measure very different oxygen atmospheres.

Le **GEN' AIR** consists of two parts:

- The pump : it depletes or enriches with oxygen the gas passing through its zirconia tube. It requires a low gas flow : between 1 and 12 l/h. It can be used for inert gas/oxygen mixtures or buffered/oxygen mixtures such as CO/CO₂/O₂ or H₂/H₂O/O₂.

- The sensor : it measures the oxygen partial pressure generated by the pump.

The use of MicroPoas®¹ gives it an excellent response time and very good measurement accuracy.

¹ - Patent ANVAR/CNRS/UNIV. Grenoble.

.ITS STRENGTHS

- Generation and analysis of atmospheres with controlled oxygen levels
- Use of carrier gases in small quantities
- Cost limitation by using a single gas
- High working dynamics
- Compact and safe system
- Limited maintenance and servicing
- Excellent measurement stability
- Measurement of oxygen partial pressure from 10⁻³⁵ to 0.25 atm



Technopole de Château-Gombert, 22-26 rue John Maynard Keynes, 13013 MARSEILLE

Tel : +33(0)4 91 95 65 12 Fax : +33(0)4 91 64 22 27 email : contact@setnag.com

www.setnag.com



.PRINCIPAL OF OPERATION



THE PUMP :

A touch-screen display allows data visualisation and pump's parameters setting. The pump can be controlled with pO₂ or current regulation.

The flow of oxygen generated by the pump in the case of oxidizing or neutral gas can be calculated using the formula :

$$X = X^0 + 0.209 \frac{I}{D}$$

Where :

X⁰ : molar fraction of oxygen before the pump

X : molar fraction of oxygen after the pump

I : current [A]

D : carrier gas flow rate [l/h]

THE GAUGE :

Placed after the pump, it is used to validate the partial pressure generated by the pump. The measure is performed using the MicroPoas®, a zirconia probe with a metallic internal reference.

Like all zirconia probes, the MicroPoas operates according to **Nernst's law** :

$$E = \frac{RT}{4F} \ln \frac{P_{mes}}{P_{ref}}$$

In the case of the MicroPoas®, the reference is set by an equilibrium a metal and its oxide.

.EXEMPLE OF PERFORMANCES

At 5l/h and 800°C, the performances obtained with different carrier gases are as follows :

Gas	pO ₂ mini	pO ₂	pO ₂ maxi
Air	17.2 %	20.9 %	25 %
Nitrogen	10 ⁻⁸ atm	10 ⁻⁷ atm	10 ⁻² atm
Ar+5%H ₂	10 ⁻³⁰ atm	10 ⁻²⁷ atm	10 ⁻²⁴ atm
CO/CO ₂ *	10 ⁻¹⁹ atm	10 ⁻¹⁷ atm	10 ⁻¹⁶ atm

* Note that the CO/CO₂ balance is not very stable in these conditions.

.TECHNICAL DATA

MEASURING PRINCIPLE MicroPoas®, zirconia probe with internal metal reference

MESUREMENT RANGE 10⁻³⁵ at 0,25 atm* O₂

FLOW RATE from 1 to 12 l/h**

OUTPUT SIGNALS RS485 link ModBus or RS232 protocol

DIMENSIONS AND WEIGHT 430x170x430 mm (wxhxp) - 15kg

POWER SUPPLY 230 Vac – 50/60 Hz

POWER 550 VA

*Measuring traces of oxygen with a zirconia probe is tricky because the presence of traces of impurities such as combustible compounds can create instability. This is particularly true in the 10⁻⁸ to 10⁻¹² atm O₂ range. The use of buffered mixtures makes it possible to generate reducing atmospheres in a controlled manner.

** Flow control must be provided by an external system. The use of a mass flow regulator is recommended (consult us).